



TI1823

Revised 4-94

Copyright, Eastman Kodak Company, 1994

## KODAK INSIGHT Thoracic Imaging Film / 4319

### 1) Description

KODAK INSIGHT Thoracic Imaging Film / 4319 is a high speed, ortho-sensitive medical x-ray film for use with special asymmetrical green emitting intensifying screens. INSIGHT Film has asymmetrical T-grain emulsion technology with absorbing dyes that reduce screen-light crossover to zero. This film features improved static protection. When used as designed with KODAK INSIGHT Thoracic Imaging Screens, KODAK INSIGHT Thoracic Imaging Film / 4319 is slightly slower (12%) than KODAK Ortho C Film / 4511 using KODAK LANEX Medium Screens, yet exhibits improved visualization of mediastinum, retrocardiac, and retrodiaphragmatic areas while yielding improved resolution of fine detail on both lung field and bony structures. It is processable in existing automated processing cycles.

This asymmetrical product is designed to be used with a specific film-screen orientation, which is facilitated by both a notch in the film and a special dye in one emulsion of the film. When the notch is in the upper right hand corner of the film, the "C" Emulsion is facing up. This emulsion must be in contact with the back, or non-tube side intensifying screen. The opposite emulsion, or "G" Emulsion, contains a special dye for darkroom identification, and must be in contact with the front, or tube side screen.

### 2) Safelight

Use a KODAK GBX-2 Safelight Filter with a frosted 15-watt bulb located at least 4 feet from the film.

Latensification: Safelight exposure after primary x-ray exposure.

Hypersensitization: Safelight exposure before primary x-ray exposure.

### 3) Storage and Handling

#### Handling -

Hands must be clean, dry and free of lotions, etc. Film should be handled carefully by the edges to avoid physical strains such as pressure, creasing, or buckling.

#### Storage -

Store unexposed film at 50 to 70 F (10 to 20 C), at 30 to 50 percent RH, and properly shielded from x-rays, gamma rays, or other penetrating radiation. Keep exposed film in a cool, dry place that is properly shielded from penetrating radiation. Process as soon as possible after exposure. Processed film should be stored at 60 to 80 F (16 to 27 C), at 30 to 50 percent RH.

#### 4) Sensitometric Parameters

Relative Speed:	Measured at a density of 1.00 above gross fog.
Lower Scale Contrast:	Measured as slope of the line between densities of 0.25 and 1.00 above gross fog.
Upper Scale Contrast:	Measured as slope of the line between densities of 1.15 and 2.50 above gross fog.
Gross Fog:	Density of film base plus processing fog.

#### 5) Process Variations

Changes to speed, contrast, and fog as a result of temperature variation from normal are included in GRAPH S Section.

#### 6) Intermix

This film can be processed with intermixes of common medical x-ray films.

Variations of bromide and iodide ions in KODAK RP X-OMAT Developer cause sensitometric speed effects that are significantly different for T-MAT Films than for conventional films; GRAPH included.

Some soiling and discoloration of the rollers of automatic processor developer racks may occur when KODAK INSIGHT Film is introduced into a film intermix.

#### 7) Automated Processing

##### Processors -

The following processors are recommended with KODAK RP X-OMAT Chemicals using the standard process cycle:

KODAK RP X-OMAT Processor, Model M7 Series

KODAK X-OMAT Processor, Model M4 Series

KODAK RP X-OMAT Processor, Model M8

KODAK X-OMAT M20 and M35 Processor

KODAK RP X-OMAT Processors, Models M5, M6

#### 8) Emergency Manual Processing

(Not recommended for regular use, but can be used when automated processor fails)

Solution/Step	Temperature	Time	Agitation
KODAK RP X-OMAT Developer working solution	80 F (26.5 C)	1 min	No agitation. Tap hanger immediately after immersion to remove film surface air bubbles.
KODAK Indicator Stop Bath OR Running Water Rinse	80 F (26.5 C)	20 sec	Continuous, moderate
KODAK RP X-OMAT X-OMAT Fixer and Replenisher	80 F (26.5 C)	1 min	Vigorous at start
Running water wash <sup>1</sup> (8 volume changes/hour)	80 F (26.5 C)	5 min	---
Dry	120 F (49 C)	---	---

<sup>1</sup> KODAK PHOTO-FLO Solution may be used after washing to minimize water spots and drying marks.

**NOTICE!** Observe precautionary information on product labels and on the Material Safety Data Sheets.

## 9) Image Structure

### Diffuse rms Granularity -

GRAPH included; read at net diffuse visual densities from 0.5 to 2.0, 48-micrometre aperture.

## 10) Graphs<sup>1</sup>

### Characteristic:

- A) RP X-OMAT Chemicals (10-91)
- B) RP X-OMAT Developer, Temperature Series (9-90)

### Reciprocity:

- C) (10-90)

### Process Variations from Normal Processing Temperature:

- D) Speed (9-90)
- E) Contrast (9-90)
- F) Fog (9-90)

### rms Granularity:

- G) (12-90)

### Safelight Sensitivity:

- H) (9-90)

<sup>1</sup>NOTICE: While the data presented are typical of production coatings, they do not represent standards that must be met by Kodak. Varying storage, exposure, and processing conditions will affect results. The company reserves the right to change and improve product characteristics at any time.

**Drying:**

I) (9-90)

**Spectral Sensitivity:**

J) (9-90)

**Bromide Effects:**

K) (9-90)

**NOTE:** The Kodak materials described in this publication for use with KODAK INSIGHT Thoracic Imaging Film / 4319 are available from dealers who supply Kodak products. You can use other materials, but you may not obtain similar results.

The contents of this publication are subject to change without notice.

**Kodak, Eastman, X-Omat, Insight, T-Mat, Photo-Flo, and Estar are trademarks.**

Health Imaging  
EASTMAN KODAK COMPANY - Rochester, NY 14650


---

**End of Data Sheet**

---

TI1823A 10-91  
CHARACTERISTIC, For Publication


KODAK InSight Thoracic Imaging Film/4319  
1/25 second Simulated Green Screen Exposure  
Seasoned KODAK RP X-OMAT Chemicals, 95 F (35 C); KODAK RP X-OMAT Processor M6;  
Diffuse Visual Densitometry



Notice: While the data presented are typical of production coatings, they do not represent standards which must be met by Eastman Kodak Company. Varying storage, exposure, and processing conditions will affect results. The company reserves the right to change and improve the product characteristics at any time.


**TI1823B 11-90**  
CHARACTERISTIC, For Publication

KODAK InSight Thoracic Imaging Film/4319  
1/25 second Simulated Green Screen;  
KODAK RP X-OMAT Chemicals KODAK RP X-OMAT Processor, Model M6;  
Diffuse Visual Densitometry



Notice: While the data presented are typical of production coatings, they do not represent standards which must be met by Eastman Kodak Company. Varying storage, exposure, and processing conditions will affect results. The company reserves the right to change and improve the product characteristics at any time.


TI1823C 11-90  
RECIPROCITY, For Publication  
KODAK InSight Thoracic Imaging Film/4319  
KODAK RP X-OMAT Chemicals, 95F,  
KODAK RP X-OMAT Processor, Model M6  
(Reference: 1/100 second = 0%)



Notice: While the data presented are typical of production coatings, they do not represent standards which must be met by Eastman Kodak Company. Varying storage, exposure, and processing conditions will affect results. The company reserves the right to change and improve the product characteristics at any time.

**TI1823D 11-90**  
TEMPERATURE VARIATION, For Publication


KODAK InSight Thoracic Imaging Film/4319  
Percent Change in Relative Speed  
KODAK RP X-OMAT Chemicals, KODAK RP X-OMAT Processor, Model M6  
(Reference: Normal Temp. = 0% Change)  
(4 F = 2.2 C)



Notice: While the data presented are typical of production coatings, they do not represent standards which must be met by Eastman Kodak Company. Varying storage, exposure and processing conditions will affect results. The company reserves the right to change and improve product characteristics at any time.

TI1823E 11-90  
TEMPERATURE VARIATION, For Publication


KODAK InSight Thoracic Imaging Film/4319  
Percent Change in Contrast  
KODAK RP X-OMAT Chemicals, KODAK RP X-OMAT Processor, Model M6  
(Reference: Normal Temp. = 0% Change)  
(4 F = 2.2C)



Notice: While the data presented are typical of production coatings, they do not represent standards which must be met by Eastman Kodak Company. Varying storage, exposure and processing conditions will affect results. The company reserves the right to change and improve product characteristics at any time.

TI1823F 11-90  
TEMPERATURE VARIATION, For Publication


KODAK InSight Thoracic Imaging Film/4319  
Density Change in Fog  
KODAK RP X-OMAT Chemicals, KODAK RP X-OMAT Processor, Model M6  
(Reference: Normal Temp. = 0% Change)  
(4 F = 2.2 C)



Notice: While the data presented are typical of production coatings, they do not represent standards which must be met by Eastman Kodak Company. Varying storage, exposure and processing conditions will affect results. The company reserves the right to change and improve the product characteristics at any time.

TI1823G 12-90  
GRANULARITY, For Publication


KODAK InSight Thoracic Imaging Film/4319  
KODAK RP X-OMAT Chemicals, 95F  
KODAK RP X-OMAT Processor, Model M6



Notice: While the data presented are typical of production coatings, they do not represent standards which must be met by Eastman Kodak Company. Varying storage, exposure and processing conditions will affect results. The company reserves the right to change and improve product characteristics at any time.

TI1823H 11-90  
SAFELIGHT SENSITIVITY, For Publication

KODAK InSight Thoracic Imaging Film/4319  
KODAK GBX-2 Safelight Filter, 15 watt lamp, 48 inches  
KODAK RP X-OMAT Processor, Model M6; KODAK RP X-OMAT Chemicals, 95 F  
(Fog growth with increasing safelight exposure)




Notice: While the data presented are typical of production coatings, they do not represent standards which must be met by Eastman Kodak Company. Varying storage, exposure and processing conditions will affect results. The company reserves the right to change and improve product characteristics at any time.

— Latensification  
- - - Hypersensitization


TI1823I 11-90  
DRYING, For Publication

KODAK InSight Thoracic Imaging Film/4319  
KODAK RP X-OMAT Processor, Model M8  
KODAK RP X-OMAT Chemicals, 96 F  
Dryer Capacity Used to Dry Film at 125 F




Notice: While the data presented are typical of production coatings, they do not represent standards which must be met by Eastman Kodak Company. Varying storage, exposure and processing conditions will affect results. The company reserves the right to change or improve product characteristics at any time.

TI1823J 11-90  
SPECTRAL SENSITIVITY, For Publication  
KODAK InSight Thoracic ImagingFilm4319  
Seasoned KODAK RP X-OMAT Chemicals; KODAK RP X-OMAT Processor, Model M6  
Effective Exposure 1.4 sec; Diffuse Visual Densitometry



Notice: While the data presented are typical of production coatings, they do not represent standards which must be met by Eastman Kodak Company. Varying storage, exposure and processing conditions will affect results. The company reserves the right to change and improve product characteristics at any time.

**TI1823K 11-90**  
BROMIDE EFFECTS, For Publication  
KODAK InSight Thoracic Imaging Film/4319  
KODAK RP X-OMAT Processor, Model M6  
Seasoned KODAK RP X-OMAT Chemicals, 95 F  
Normal Level is 3.5 g/L



Notice: While the data presented are typical of production coatings, they do not represent standards which must be met by Eastman Kodak Company. Varying storage, exposure and processing conditions will affect results. The company reserves the right to change and improve product characteristics at any time.